

Aligner consumables

Packaging, labeling & accessories



Zip lock bags



Label printer



Elastics



Tray seaters



Cleaning tablets



Aligner scissors

Packaging & labeling

A minimal set to neatly assemble and label an aligner kit

01

Zip lock bags

Zip Lock Bags

Size: 10.2 × 12.7 cm
or 10 × 15 cm

**02**

Label printer

Dymo LabelWriter 550

For fast labeling

**03**

Labels / tape

89 × 36 mm

Standard size for
aligner bags

**04**

Aligner case

dental case /
retainer case

For storage and
patient delivery

**05**

Elastics

Elastics Rings 3.5oz

Most often 1/8 and
3/16 inch

**06**

IPR gauge

IPR Gauge

For interproximal
reduction



For patient and clinician

Accessories for cleaning, seating, aligner removal, and precise clinical work

07

Cleaning tablets

retainer cleaning tablets

For aligner hygiene



08

Tray seaters

tray seaters

Help seat aligners fully



09

Removal hook

Aligner hook

Easy aligner removal



10

Orthodontic hook

orthodontic hook

For additional clinical tasks



11

Elastic hooks

Orthodontic elastics hook

For removing and placing elastics



12

Aligner scissors

scissor surgical curved

Choose curved scissors, 14 cm or longer

